


la belle CL 210



la belle CL 210




- BTE-module for eyeglasses, cold bendable
- digital hearing system
- natural sound transfer via ex-receiver system *
- 3 programmes selectable by push button
- 2 channel WDRC
- Dynamic Contrast Detection
- 12-band gain adjustment
- MPO
- adaptive feedback cancellation (until -17dB)
- programme tones and battery tones
- ideal for open-ear-fitting
- battery type 312
- replaceable cerumen defenser filter
- tube warm bendable
- colours: gold-brown, silver-black, anthracite-black
- mono, stereo, CROS, BiCROS

*ex-receiver = receiver in the ear canal



Technical data

		DIN EN 60118-0	DIN EN 60118-7
Max. output level OSPL90	(db)	119	110
Max. gain 50dB(FOG)	(dB)	49	42
Frequency range	(Hz)	200-6000	200-6500
Equivalent input noise	(dB)	22	22
Operating current	(mA)	1,2	1,2



Before programming, please remove the device cover plate.

Programming:

Bruckhoff fitting software and Siemens programming cable (CONNEXX cable grey)

Technical change without prior notice