

freestyle CL 210 + tinnitus noiser + volume control




Hearing system with noise generator for tinnitus therapy as one device.

Adjustable white band noise to your own personal needs with volume control (VC)




freestyle CL 210 with option TID

- modern digital communication system
 - adjustable noise level by volume control (VC)
 - 3 programmes selected by push button
 - 2-channel WDRC
 - Dynamic Contrast Detection
 - 12-band gain adjustment
 - MPO
 - adaptive feedback cancellation (until -17dB)
 - programme tones and battery tones
 - battery type 312
 - tube warm bendable
 - replaceable cerumen defender filter
 - can be adjusted to fit any individual ear shape
 - natural sound transfer quality via ex-speaker system*
 - colours: titanium/transparent
- * ex-speaker = speaker direct in the ear canal



Technical data

		DIN EN 60118-0	DIN EN 60118-7
Max. output level OSPL90	(dB)	119	110
Max. gain 50dB(FOG)	(dB)	49	42
Equivalent input noise	(dB)	22	22
Frequency range	(Hz)	200-6000	250-6500
Operating current	(mA)	1,2	1,2



Programming:
with bruckhoff fitting software and Siemens programming cable (CONNEXX cable grey)

Technical change without prior notice. Before programming, please remove the device cover plate